



## 13

## LIPID STAIN

## 13.1 INTRODUCTION

The Oil Red O (ORO) stain can identify neutral lipids and fatty acids in smears and tissues. Fresh smears or cryostat sections of tissue are necessary because fixatives containing alcohols, or routine tissue processing with clearing, will remove lipids. The ORO is a rapid and simple stain.



## OBJECTIVES

After reading this lesson, you will be able to:

- explain the principle of lipid stain
- describe various reagents used for lipid stains
- describe the procedure of lipid staining.

## 13.2 LIPID STAIN

**Aim:** To demonstrate intracellular lipid in tissue sections.

**Principle:** The dye being more soluble in the lipid to be demonstrated than in the vehicular solvent. The polyazo group of dyes includes the oil red series, the sudan red series and sudan blacks. All these dyes are interchangeable and may be substituted.

Sudan series - Sudan III  
- Sudan IV  
- Sudan black

**Control** - Lipid positive section



## Lipid Stain

### Note

- Use cryostat sections of 8 to 10 micron thickness or formalin fixed smears.
- working Oil Red O solution should be freshly prepared from stock solution and kept in close container.
- Never take the sections through clearing solvent prior to mounting as this will remove the lipid to be demonstrated.
- Frozen sections should be used to stain neutral triglycerides.
- Lipoproteins may be demonstrated on paraffin sections.
- Alcohol fixation removes most lipids.



### INTEXT QUESTIONS 13.1

1. Dyes used in lipid stain is .....
2. .... is used to stain nuclei
3. Lipid is demonstrated by ..... colour
4. Nuclei is demonstrated by ..... colour
5. Lipoproteins may be demonstrated on ..... section
6. .... fixation removes lipids



### WHAT HAVE YOU LEARNT

- Lipid stain is used to demonstrate intracellular lipid in tissue sections
- Polyazo group of dyes like oil red series, sudan red series and sudan blacks are the dye used for demonstrating lipids in tissue section
- The frozen section should be cut 8 to 10 micron thickness
- Lipids appear as red and nuclei appear blue
- Lipoproteins may be demonstrated on paraffin sections
- Alcohol fixation removes most lipids

## MODULE

Histology and Cytology



Notes

## MODULE

Histology and Cytology



Notes

Lipid Stain



### TERMINAL QUESTIONS

1. What is principle of lipid stain?
2. Name three dyes used to demonstrate lipid in tissue sections.
3. What precautions should be observed during lipid staining on tissue sections?
4. What should be the thickness of sections for lipid staining?
5. What is the mounting media used in lipid staining?



### ANSWERS TO INTEXT QUESTIONS

#### 13.1

1. Sudan series
2. Alum hematoxylin
3. Red
4. Blue
5. Paraffin
6. Alcohol